

QA-1290 Plug-In Read Me

FLUKE.

Biomedical

ansur QA-1290 Plug-In is a product facilitating automatic test and inspection with the use of the QA-1290 NIBP-Analyzer.

Introduction

This read me file contains last minute information for Fluke Biomedical's ansur QA-1290 version 2.0.0.

History

This chapter describes the changes made between each new released version.

Release 2.0.0

- Altered names and graphics for Fluke version.
- Added support for ansur V2.

Release 1.2.7

- Added French language support
- Fixed losing communication during test when dragging the Test Guide.
- Fixed bug pressing next and previous, this no longer demands user to re-execute test
- Fixed setting from 200/14 to 200/140
- Fixed Enter key behavior when pressing Enter while running a test
- Added Readme file to the QA-1290 installation set Readme(QA1290).htm

Fluke Corporation

PO Box 9090
Everett, WA USA 98206

Fluke Eurpoe B.V.
PO Box 1186, 5602 BD
Eindhoven, The Netherlands

For more information call:

U.S.A. (800) 443-5853 or Fax (425) 356-5116
Europe / M-East (31 40) 2 678 200 or Fax (31 40) 2 678 222
Canada (905) 890-7600 or Fax (905) 890-6866
Other countries (425) 356-5500 or Fax (425) 356-5116
Internet: <http://www.fluke.com>

© 2005 Fluke Corporation. All rights reserved.